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Introduction

Flu is a potentially serious viral respiratory disease caused by influenza virus that circles the globe every year. As it spreads, the virus slowly changes its genetic materials that make it very difficult for our immune systems to recognize the infection. Three types of Influenza virus have been identified and are named influenza A, B and C that are divided into subtypes based on two glycoproteins on the surface of the virus, the haemagglutinin (HA) and neuraminidase (NA). There are fifteen different HA (i.e. H1-H15) and nine different NA (i.e. N1-N9). The current subtypes of influenza A viruses that infect humans are new H1N1. The influenza virus originally came to humans from birds. Recently, a large outbreak of novel influenza A virus, is originated from swine, jumped the species barrier and began to infect people.

These rapid influenza diagnostic test kit can identify Influenza A&B virus directly from nasal swab or Nasopharyngeal swab with a high degree of accuracy.

Material provided Influenza A&B kit contains following items

- 1) Influenza A & B test strip wrapped with foil pouch
- 2) Sterilized swab for specimen collecting
- 3) Buffer for specimen extraction
- 4) Instruction for use

Indications

For qualitative detection Influenza A & B virus antigen from human nasal, nasopharyngeal and throat swab.

Specimen collection (refer to image 1)

1. Because collecting sample has a chance of infection unknown, it should be treated by trained person having special technique.
2. As collecting sample affects to test result crucially, it is needed to be careful.
3. When sample collection is done at the several part(Nasal and Throat), can get more exact result.
4. Specimen should be tested as soon as possible after collection.

< 1 > Nasal swab specimens

To collect a nasal swab specimen, it is very important to collect sample from both nostrils. Insert the inclosed swab into a nostril with care. Rotating softly, push the swab until resistance is met at the level of turbinate. Collect the exudate on the nasal wall. And then collect the exudate from the other nostril with the swab just used at the same way.

< 2 > Nasopharyngeal swab specimens


To collect nasopharyngeal specimen, insert a sterilized swab into a nasal cavity deeply to touch nasopharyngeal. Rotating softly, follow the instruction on nasal swab specimen collection.

< 3 > Throat swab specimens

To collect a throat swab specimen, vigorously rub a swab on both tonsillar surfaces and the posterior pharynx.

Specimen Collection

(Image 1 Specimen collection)

Storage of collected swab and sample


When the sample have to be stored, swab sample itself could be stored in the cold storage (2~8℃) for 12hours.

Test procedure (refer to image 2)

1. This test should be done at the room temperature, so when test reagent was stored in the cold storage, before doing the test 15~30minutes, place the reagent in the room temperature and make it same as room temperature. When reagent stored at the room temperature, can open and use it immediately.
2. Insert the patient swab sample into the test tube, then, swirl the swab 5~10 times while pressing the head against the bottom and side of the test tube.
3. Roll the swab head against the inside of the tube while you removing it. Dispose the used swab in accordance with your biohazard waste disposal protocol.
4. Take the test strip from the foil pouch.
5. Place the test strip into the extraction tube as arrow indications on the test strip pointing down. Do not handle or move the test strip until the test is completed and ready for reading.
6. Interpret the result between 15~20 minutes.
7. Please do not interpret over 20 minutes.

Test Procedure

(Image 2 Test Procedure)

Interpretation of the test (refer to image 3)

1. Negative: Only one band appears on control line(C).
2. Influenza A Positive(A) : one purple line(C), one gray line(A)
A means gray line which is closely located in the bottom.
3. Influenza B Positive(B) : one purple line(C), one gray line(B)
B means gray line is closely located in the middle.
4. Influenza A & B Positive : one purple line(C), two gray line(A,B).
5. Invalid result : If the control band fails to appear within the result window, the result is considered invalid. The directions may not have been followed correctly or the test may have deteriorated. It is commended that the specimen be retested.

Interpretation

(Image 3 Interpretation of the result)

Performance characteristics

1. Sensitivity and Specificity

Nasal Swab

		A	B	Neg.	Total
P C R	A	55	0	3	58
	B	0	26	2	28
	Neg.	0	0	101	101
	Total	55	26	106	187

A type sensitivity : 94.82%

B type sensitivity : 92.85%

Specificity : 100%

Nasopharyngeal Swab

		A	B	Neg.	Total
P C R	A	52	0	3	55
	B	0	27	3	30
	Neg.	0	0	80	80
	Total	52	27	86	165

A type sensitivity : 94.54%

B type sensitivity : 90%

Specificity : 100%

2. Cross reactivity

- Type A influenza virus

Cross-reactivity at the B test line position was not observed with the following type A influenza viruses (10^5 to 10^8 pfu/ml)

A/New Caledonia/20/99(H1N1), A/California/7/2009 (H1N1), A/Panama/2007/99(H3N2), Influenza A/Anhui/1/05 (H5N1)

- Type B influenza virus

Cross-reactivity at the A test line position was not observed with the following type B influenza viruses (10^5 to 10^8 pfu/ml)

B/Yamanashi/166/98, B/Panama/45/90, B/Brisbane/60/08, B/Massachusetts/02/2012, B/Yamagata/16/88, B/Phuket/3073/2013, B/Malaysia/2506/2004, B/Wisconsin/1/2010, B/Victoria/2/87, B/Tokio/53/99, B/Victoria/504/00, B/Malaysia/2506/04, B/Florida/07/04, Florida/04/06

- The other viruses

Cross-reactivity was not observed with any of the following viruses (10^5 to 10^8 TCID₅₀/ml)

Adenovirus (type 1,2,3,4,5,6,7), Parainfluenza virus (type 1,2,3,4), Respiratory syncytial virus (type A, B), Rhinovirus (type A, B)

Limitation of the test

This kit is developed for detection of Influenza type A & B virus. This kit can provide rapidity and easiness for the test. But it still remains the possibility of false positive or false negative result caused by various factors. So other clinically useful tests are required for accurate test. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test. The results must be evaluated in conjunction with other clinical data available to the physician.

Warnings

1. For *in vitro* diagnostic use only.
2. Do not RE-USE all materials.
3. Wear protective gloves while handling specimens. Wash hands after doing test.
4. The test device is sensitive to humidity and heat. Perform the test immediately after removing the device from the foil pouch.
5. Do not use the device beyond the expiration.
6. Decontaminate and dispose of all specimens, tested devices and potentially contaminated materials, as if they were infectious waste, in a biohazard container.

Package

25 tests / 1 Kit

Kit storage and stability

The test kit should be stored at room temperature.

This product is stable for 24 month since the date of manufacture.

Manufactured by - RapiGEN INC
2F, 25, Heungan-daero, Gunpo-si, Gyeonggi-do, Republic of Korea
Tel : +82-31-427-4677, FAX : +82-31-427-4678, Web : www.rapigen-inc.com
Altenhofstrasse 80, D-66386 St.Ingbert, Germany
Tel : +49-6894-581020, FAX : +49-6894-581021, Email : info@mt-procons.com www.mt-procons.com

G10RHS25-K/E04(2017.10.10)